VIAN

Vision beyond limits

Vian Biotech Private Limited



Eyes are precious
Deserve to be CARED

Scanned with OKEN Scanner

PRODUCTS GROUP	PAGE NO.
COMPANY PROFILE	2
IOLs	3-6
EYE DROPS	7-8
B-HEX PUPIL EXPANDER & FORCEPS	9
OPHTHALMIC PRODUCTS & CONSUMABLES	10
INSTRUMENT CARE	11
BIMANUAL SETS & HANDLE	12
BIPOLAR CABLE & FORCEPS	13
BLADE BREAKERS & NEEDLE HOLDERS	14 - 16
CANNULAES HANDLE	17 - 20
CAUTERYS, LOOPS & SPATULAS	21 - 22
CHOPPERS	23
CLAMPS	24 - 25
CURETTES, PROBES, SPOONS & DILATORS	26 - 27
DISSECTORS & CLAMPS	28
FORCEPS	29 - 35
HOOKS & POLISHERS	36
HOOKS	37 - 38
I.O.L. MANIPULATORS	39 - 40
KERATOMES KNIVES & SPUDS	41 - 42
MALLETS & FORCEPS	43
MARKERS FIXATION RING	44
PUNCHES & TREPHINES	45
RHEXIS, HOLDER & FOLDER	46 - 48
RETRACTORS & KNIVES	49
SCISSORS	50 - 54
SPECULUMS	55 - 57
SPATULAS, CALIPERS & MARKERS	58
STERILIZATION CASE	59
VITREO-RETINAL INSTRUMENT & PEELARS	60

Scanned with OKEN Scanner

COMPANY PROFILE

VIAN was started in 2022 with vast experience in the field of ophthalmology sales & marketing of medicines / micro surgical instruments for more than 10 years.

VIAN committed to provide high quality Titanium and Stainless steel surgical instruments and trading of applied ophthalmic items with the reasonable team of experienced sales, marketing and administrative staff with intent to provide the qualitative & innovative formulation to help protect, preserve and restore the vision in patents through predictable outcomes and improved compliance.




VISION

To become most trusted and preferred brand for our partners, HCPs and patients.

MISSION

Grow together, by involving Healthcare professionals and partners to provide innovative, Quality products and services to patient and communities.





VALUES

Vital: To be Vital, full of energy and lively
Innovative and unique in the solutions and services we provide with integrity
Aspired to fulfill 100% satisfaction on commitment what we give.
Nascent - To be a future potential along with out partners and people whom we involve.

Scanned with OKEN Scanner

VIAN
Vision beyond limits

ifocus® IOLs

PHOBIC

HA 6030

Single piece Sterile with 360° Square Edge **HYDROPHOBIC ASPHERIC** Acrylic Foldable Aspheric Posterior Chamber Intraocular Lens (IOL) with UV blocking, Optionally with blue light filter.

This optic is **MONOFOCAL**

The IOL is non-preloaded, which has to be manually loaded into a compatible Injector.



Hydrophobic Acrylic Foldable - IOL

POWER RANGE	
Standard Powers	+1.0 D→+30.0 D (0.5 D steps)

MATERIAL

Hydrophobic Acrylic blanks contain an UV Absorber with <0.5% Water Content.

Refractive Index	1.56
Sterilization	Ethylene Oxide

DESIGN

Optic Diameter Overall Diameter	6mm	13.00mm
Optic Design	Aspheric-Aberration Neutral Biconvex (+10.0D to 35.00 D)	
Nominal A-constant	118.07	
Haptic Design	Modified 'C' Haptic	
PCO Prevention	360° Square Edge	
Shelf Life	3 years after sterilization	

INJECTION

Estimated Incision Size	2.40mm / 2.6mm
Injector & Cartridge	Glojet - Click lock Mechnism

CONSTANTS
for optical and immersion ultrasound measurements

SRK/T A-Constant	SRK/II A-Constant
118.4	118.7

Scanned with OKEN Scanner

VIAN
Vision beyond limits

ifocus® IOLs







	fold Hydrophilic spheric Foldable-IOL	AO fold Hydrophilic Aspheric Foldable-IOL	BLF fold Hydrophilic yellow Aspheric Foldable-IOL
Standard Powers	00.0 D→+30.0 D (0.5 D steps)	00.0 D→+30.0 D (0.5 D steps)	00.0 D→+30.0 D (0.5 D steps)
Extreme Powers	-10.0 D→-1.0 D (0.5 D steps) +31.0 D→+45.0 D (0.5 D steps)	-10.0 D→-1.0 D (0.5 D steps) +31.0 D→+45.0 D (0.5 D steps)	-10.0 D→-1.0 D (0.5 D steps) +31.0 D→+45.0 D (0.5 D steps)
MATERIAL	Copolymer of hydroxyethyl methacrylate and methyl methacrylate with 26% water content, UV blocker	Copolymer of hydroxyethyl methacrylate and methyl methacrylate with 26% water content, UV blocker	Copolymer of hydroxyethyl methacrylate and methyl methacrylate with 26% water content, UV blocker with blue light filter
Refractive Index	1.56	1.56	1.46
Sterilization	Steam	Steam	Steam
Nominal A-constant	118.0	118.0	118.0
Haptic Design	0° Angulation	0° Angulation with Posterior Vaulting	0° Angulation with Posterior Vaulting
PCO Prevention	360° Square Edge	360° Square Edge	360° Square Edge
Shelf Life	3 years after sterilization	3 years after sterilization	3 years after sterilization
Estimated Incision Size	2.20mm	2.00mm to 2.80mm	2.00mm to 2.80mm
Injector & Cartridge	Glojet - Click lock Mechnism	Glojet - Click lock Mechnism	Glojet - Click lock Mechnism
CONSTANTS	*for optical and immersion ultrasound measurements*	*for optical and immersion ultrasound measurements*	*for optical and immersion ultrasound measurements*
A-Constant	A = 118.0*	A = 118.0*	A = 118.0*

Scanned with OKEN Scanner


VIAN
Vision beyond limits

ifocus®
Yellow PMMA Intraocular Lens

Yellow PMMA Aspheric Optics with 360° Square Edge Intraocular Lens

GPY 6025


Standard Powers	00.0 D→+35.0 D (0.5 D steps)
Extreme Powers	-10.0 D→-1.0 D (1.0 D steps) +35.0 D→+45.0 D (1.0 D steps)
Optic Diameter **Overall Diameter**	6.00 mm 12.50 mm
Optic Design	Aspheric with Square Edge,
Haptic Design	Modified "C" StepVault
Est Ac Depth	5.00 mm
Dialing Holes	2
Suggested A Constant	118.4*
Shelf Life	5 years after sterilization



ifocus®
IRIS CLAW Intraocular Lens

PS 5580

Standard Powers	+2.00 to +30.0 D
Optic Diameter **Overall Diameter**	5.50 mm 8.00 mm
Optic Design	Bi-convex
Ac Depth	0 mm
Suggested A Constant	Anterior Side - 115.0* Posterior Side - 117.2*
Shelf Life	5 years after sterilization



Scanned with OKEN Scanner

VIAN
Vision beyond limits

ifocus® PMMA IOLs

	UPS 5220 Posterior Chamber - IOL	UPS 6025 Posterior Chamber - IOL	UPS 6025 Anterior Chamber-IOL	UPS 6535 Sclearal Fixction-IOL
Standard Powers	00.0 D→+35.0 D (0.5 D steps)	-10.0 D→+35.0 D (0.5 D steps)	00.0 D→+35.0 D (0.5 D steps)	10.0 D→+35.0 D (0.5 D steps)
Extreme Powers	10.0 D→1.0 D (1.0 D steps)	10.0 D→1.0 D (1.0 D steps) +30.0 D→+40.0 D (1.0 D steps)	-10.0 D→+40.0 D	-10.0 D→+40.0 D
Optic Diameter **Overall Diameter**	5.25 mm 12.00 mm	6.00 mm 12.50 mm	6.00 mm 12.50 mm	6.00 mm 12.50 mm
Optic Design	Equiconvex	Equiconvex	Equiconvex	Equiconvex
Haptic Design	"C" with StepVault	"C" with StepVault	"S" with StepVault	"C" with StepVault
Ac Depth	5.00 mm	5.00 mm	2.95 mm	4.8 mm
Dialing Holes	Nil	2	Nil	2
Suggested A Constant	118.4*	118.4*	115.3*	118.2*

MED

3 Piece
Hydrophobic & PMMA 3 Piece-IOL

Standard Powers	10.00 D → +30.0 D (0.5 D steps)	**Optic Diameter** **Overall Diameter**	6 mm 13.00 mm, 13.50 mm
Haptic Material	PVDF (Polyvinylidene	**Optic & Haptic Design**	Biconvex - 360° Square Edge with Aspheric Optics & Modi - 0° Angulation
Refractive Index	1.493	**Nominal A-constant**	118.8
Sterilization	Ethylene Oxide	**PCO Prevention**	360° Square Edge
Estimated Incision Size	2.6 - 2.8 mm	**SRK/T A-Constant** **SRK/ II A-Constant**	A = 118.7* & A = 118.9*
Injector & Cartridge	Med Jet™-Front Load with 2.8mm Cartridge		
MATERIAL	Hydrophobic Acrylic blanks contain an UV Absorber with <0.5% Water Content. optional blue light tering Chromophore		

Scanned with OKEN Scanner

VIAN
Vision beyond limits

EYE DROPS

Look Fresh™ EYE DROPS Tears

- Provides relief from burning irritation & discomfort.
- Helps to soothe the irritation and burning seen in dry eyes by maintaining proper lubrication of the eyes.
- Has other osmoprotectants Glycerin, Erythritol, Levocatnitine, Calcium Chloride, Magnesium Chloride, Potassium Chloride.
- Osmoprotectants can reduce corneal staining, decrease cell apoptosis and inflammatory cytokines and increase the number of goblet cells.

Hylovian™ Eye Drops

- These may be caused by your eye not producing enough tears to keep the eye moist.
- Provides soothing relief from the symptoms of dry eyes, such as soreness, irritation or a gritty sensation.
- It is Preservative-free and can be used by people who wear rigid or gas permeable contact lenses, or who are sensitive to preservatives
- Sodium hyaluronate has also been shown to stimulate the development of epithelial cells on the eye's surface, which can help to heal damage to the surface of the eye resulting from lack of lubrication.

Look Fresh™ EYE DROPS Gel

- Excellent lubricating efficacy, keeps eyes cool & moist
- Superior electrolytes balance increases the goblet cell density
- Superior mucoadhesive force of 192.4 prolongs the drug retention time
- Has other osmoprotectants Glycerin, Levocarnitine, Erythritol, Calcium Chloride, Magnesium Chloride, Potassium Chloride.

Pegafresh™ EYE DROPS

- Provides long lasting relief form dry eyes and associated burning, itching, irritation and discomfort
- Controls and stabilizes deficient tear production
- Thickens the tear film and prolongs tear film break up time.
- Produces better best-corrected visual activity scores as compared to CMC/Glycerin
- Has Hp Guar Gum, Sorbitol, AMP Ultra PC 2000, Boric Acid, Potassium Chloride, Sodium Chloride, Edetate Disodium.

Predvian™ EYE DROPS

- It is a broad-spectrum antibiotic that is active against both Gram + and Gram- bacteria.
- used only to treat / prevent infections that proven or strongly suspected / caused by bacteria.
- Prednisolone ophthalmic solution is indicated for steroid-responsibe inflammatory ocular conditions for which a corticosteroid is indicated and where bacterial infection or a risk of bacterial ocular infection exists.
- The combination can also be used for post-operative inflammation and any other ocular inflammation associated with infection.



Scanned with OKEN Scanner

EYE DROPS

VIAN
Vision beyond limits

Moxivian™ EYE DROPS

- Broad sprectrum coverage
- Improved Ocular Bioavailability
- Reduces probability for developing bacterial resistance
- Reducing ability of cell's efflux pump to flush out antibiotic
- Increases drug's stay in bacterial cells for longer duration of action.

Moxivian™-DX Eye Drops

The Fixed Combination of Moxifloxacin & Dexamethasone was found to be well tolerate & effective in minimizing inflammation following Cataract Surgery

Moxivian™-LP Eye Drops

GREATER SAFETY MARGIN

- Seasonal Allergic Conjunctivitis
- Low risk of IOP elevation
- Low risk of steroid induced cataract formation
- Excellent choice for long term maintenance therapy

GREATER EFFICACY

- Higher Lipophilicity for greater penetration

ATROVIAN EYE DROPS एट्रोवियन

Key Uses of Atropine Eye Drops:

- **Eye Examinations:** Dilates the pupil to allow easier viewing of the back of the eye.
- **Inflammation (Uveitis/Iritis):** Reduces pain and inflammation by relaxing muscles, helping the inflamed tissue recover.
- **Amblyopia (Lazy Eye):** Used in children to help strengthen the affected eye.
- **Myopia Control :** Slows the progression of nearsightedness in children and younger teens.
- **Cycloplegia:** Temporarily paralyzes the focusing muscle for diagnostic purposes, such as checking for blurred vision.

Nepavian™ EYE DROPS

NEPAVIANZ-OD EYE DROPS

- Onset of action is within 15 minutes and persist for more than 8 hours
- Corneal permeability Co-efficient of Nepafenace is 4, 19 and 28 times greater than Diclofenac, Bromfenac & Ketorolac respectively
- Effective in treating pain and inflammation associated with cataract surgery
- Effective in treatment of cystoid macular edema
- Safest Preservative SOC




Scanned with OKEN Scanner


Vian Biotech Pvt Ltd

US FDA Registered


B-HEX® Pupil Expander



B-HEX® Pupil Expander 2.0



B-OCTA™ 8.00 Pupil



Pupil Expansion Redefined


The **'B-HEX Pupil Expander'** (6.50mm diameter, 0.075 mm thick), **'B-HEX Pupil Expander 2.0'** (6.75 mm diameter, 0.125mm thick), & **'B-OCTA 8.00 Pupil Expander'** (8.00mm diameter, 0.125mm thick), are US FDA Regd. & patented products from Med Invent Devices Pvt. Ltd. These are sterile single use flexible plastic devices with notches and flanges disposed in a single plane. Alternate flanges are tucked under the iris to engage the notches to the pupil margin and provide an expanded pupil. The **B-HEX** & **B-HEX 2.0** are used only for **Phacoemulsification (Phaco) surgery.** The **B-OCTA 8.00** is used for **Manual Small Incision Cataract Surgery (MSICS) & Phaco surgery.** Each device is loaded in a housing with a transparent lid which delivers the device sterile at the incision. Each unit is placed individually in a sterile pack. The device is inserted, engaged to the pupil margin, and removed using a **'B-HEX 23 G forceps'.**

B-HEX Videos

Scan or Click

B-OCTA Videos

B-HEX® 23 G Forceps



Regd. Trademarks

B-HEX® & MED INVENT DEVICES®

The **'B-HEX 23 G forceps'** is a US FDA registered, ergonomically designed, non-sterile, reusable, Titanium **Microsurgical forceps** with a short curved stalk and serrated jaws. The jaws open in a plane perpendicular to the squeeze handle. It is a necessary accessory used for grasping and manipulating the 'B-HEX & B-OCTA Pupil Expanders'. It is sterilized by Moist heat (Autoclaving).

Authorized Distributor

Vian Biotech Pvt Ltd

534, Water Tank Road, Gattigere, Rajarajeshwari Nagar, **Bengaluru** - 560098.

Ph: **97405 26706 / 99168 60403**, info@vianbiotech.com

Scanned with OKEN Scanner

OTHER **OPHTHALMIC** PRODUCTS

BYNOCS - Amblyopia Management / Treatment Program

- Aniridia Implant (Ophtec)
- CT Rings
- Iris Hooks
- Capsule Hooks
- Duvet Wet wipes
- I Wipes
- Welch Allyn Battery
- Artificial Eye Box

- Iris Retractor of 5
- Trephins
- Artificial Anterior
- Tonometer
- Endothelial Punch
- Eye Ball Stand

- Eye Patches
- Opthalmic Strips
- Trail Lens set
- Toric Marker set
- Lower Vision Kit
- Vision Therapy Kit

iCredance SURGICAL **CONSUMABLES**

- Cresent Blade
- Keratome 2.8 Blade
- Keratome 3.2 Blade
- Side port Blade
- MVR Blades
- Surgical Drapes (All Sizes available)
- Trolly Sheets
- Black & White Goggles
- Surgeon gown
- Reading glasses
- Surgical sutures
- Eye Models



Scanned with OKEN Scanner

VIAN
Vision beyond limits

INSTRUMENT CARE

CLEANING

- Distilled or mineral water (Lab tested soft water) should be used to Decontaminate & rinse the instruments.
- For removal of contaminated stuck on the surface, use soft bristled brush or ultrasonic cleaner.
- Cannulas and hand pieces should be cleaned with distilled water.
- Then use syringe or compressed dry air to remove wet.

STERILIZATION


- Use distilled water in autoclave.
- Autoclave should be rustproof from inside otherwise rust may transfer by steam.
- During dry heat sterilization the temperature should be fully controlled otherwise instruments hardness may get affected.

STORAGE


- Instrument should be fully dried by hair dryer or in oven with temperature at 70-80OC. Medical grade oil can be used on serration joints locks, hinges & Knurled parts of the instruments.

Scanned with OKEN Scanner


BIMANUAL SETS AND HANDLE



DR. PARUL SHAH'S
Type bimanual set
VB-2202
Simcoe cortex 22gauge
Teplon handle knob
With tubing



DR. J.K. SHAH'S
Type bimanual set
VB-2208
Teplon handle
Simcoe cortex
1/A (23g) wash
Cannula (19g)



I/A SYSTEM
Irrigating aspirating
ss with titanium tip
VB-2214

BIMANUAL
Handpiece for bimanual
technique
VB-2220 (A) Irrigating
VB-2220 (B) Aspirating



B.P HANDLE
VB-2226

Scanned with OKEN Scanner

BIPOLAR CABLE & FORCEPS




Bipolar Cable
VB-4028




Bipolar Forceps
VB-4034

Scanned with OKEN Scanner


BLADE BREAKERS & NEEDLE HOLDERS




SWISS
Model blade holder
pencil 12 cm.
VB-1722



BARRAQUER
Blade breaker & holder
VB-1734



KALT
Needle holder
VB-1740 (A) Straight
VB-1740 Curved



SILCOCK
Needle holder
VB-1746

Scanned with OKEN Scanner

BLADE BREAKERS & NEEDLE HOLDERS



CASTROVIEJO

Needle holder standard jaws
VB-1752 (A) Straight Without Lock
VB-1752 (B) Straight With Lock
VB-1752 (C) Curved Without Lock
VB-1752 (D) Curved With Lock



Mc PHERSON

Needle holder delicate jaws
VB-1770 (A) Straight Without Lock
VB-1770 (B) Straight With Lock
VB-1770 (C) Curved Without Lock
VB-1770 (D) Curved With Lock



BARRAQUER

Needle holder standard model
VB-1776 (A) Without Lock
VB-1776 (B) With Lock



BARRAQUER

Needle holder
VB-1782 (A) Without Lock
VB-1782 (B) With Lock



BARRAQUER

Micro needle holder medium
VB-1794 (A) Without Lock
VB-1794 (B) With Lock

Scanned with OKEN Scanner

BLADE BREAKERS & NEEDLE HOLDERS



BARRAQUER
Micro needle holder short model
VB-1800 (A) Without Lock
VB-1800 (B) With Lock



BARRAQUER
Micro needle holder extra delicate
VB-1806 (A) Straight Without Lock
VB-1806 (B) Straight Lock
VB-1806 (C) Curved Without Lock
VB-1806 (D) Curved Lock



NEEDLE HOLDER
Micro extra delicate
VB-1812 (A) Straight Without Lock
VB-1812 (B) Straight Lock
VB-1812 (C) Curved Without Lock
VB-1812 (D) Curved Lock



TROUTMAN
Micro needle holder
VB-1818 (A) Without Lock
VB-1818 (B) With Lock



MCINTYRE
Needle holder
VB-1824 (A) Without Lock
VB-1824 (B) With Lock



LIMS
Needle holder
VB-1830 (A) Without Lock
VB-1830 (B) With Lock

Scanned with OKEN Scanner

CANNULAES HANDLE



LACRIMAL
Cannula, 23gauge
VB-1836 (A) Straight
VB-1836 (B) Curved



CANNULA
A.C. Wash curved 20gauge
VB-1842



BISHOP-HARMON
A.C. Wash cannula 10mm
19 gauge
VB-1848



RANDOLPH
Cyclodialysis cannula
curved spatulated tip
angled at 13mm,
VB-1854 (A) 19gauge
VB-1854 (B) 21gauge



CASTROVIEJO
Cyclodialysis cannula
top opening
VB-1860 (A) 10mm
VB-1860 (B) 15mm



KELMAN
Cyclodialysis cannula
spatulated tip side
opening 19gauge
VB-1866 (A) right
VB-1866 (B) left



SALZ
Aspirating cannula
15mm long 23gauge
straight thin wall
VB-1872



GILLS
Aspirating cannula
smooth blunt tip
angled 25 gauge
VB-1878 5mm tip
VB-1884 7mm tip
VB-1890 8mm tip
VB-1896 10mm tip



MORRISON
Aspirating cannula
smooth blunt tip
angled at 7mm 23gauge
VB-1902

VISCOELASTIC
Aspirating cannula
angled at 10mm from tip
VB-1908 (A) 19 gauge
VB-1908 (B) 22 gauge



BRACKEN
Irrigating cannula
smooth beveled tip
front opening 19 gauge
VB-1914



BINKHORST
Irrigating cannula
smooth blunt tip
front opening 22 gauge
VB-1920 (A) Straight
VB-1920 (B) Right
VB-1920 (C) Left



WILBRANDT-BINKHORST
Irrigating cannula smooth
blunt tip front
opening 23 gauge
VB-1926 (A) Straight
VB-1926 (B) Right
VB-1926 (C) Left



MCINTRYE-BINKHORST
Irrigating cannula smooth
blunt tip front
opening 26 gauge
VB-1932 (A) Straight
VB-1932 (B) Right
VB-1932 (C) Left



WEST
Irrigating cannula
blunt tip, 23 guage
side opening
VB-1938

Scanned with OKEN Scanner

CANNULAES HANDLE



KNOLLE-PEARCE
Irrigating vectis three
front opening ports
23 gauge loop
VB-1944



'J' SHAPED
Irrigating aspirating Cannula
VB-1986 (A) 12 'o' clock
VB-1986 (B) Right
VB-1986 (C) Left



DREWS
Irrigating vectis
one front and two
backward opening
parts, 23 gauge loop
VB-1950



Mc INTYRE
Inner cannula 23gauge
VB-1992



SIMCOE
Irrigating lens
nucleus loop two front opening
ports finely, serrated, strong
curved 25 gauge
VB-1956 (A) Right
VB-1956 (A) Left



Mc INTYRE
Outer cannula straight
blunt tip front
opening 19gauge
VB-1998




SHEETS
Irrigating vectis three
large front opening ports
serrated 20 gauge loop
VB-1962



Mc INTYRE
Anterior chamber cannula
smooth blunt tip
VB-2004 23gauge angled
VB-2010 26gauge angled
VB-2016 30gauge angled



SIMCOE
Irrigating aspirating
cannula
with silicone tubing
VB-1968A-21G
VB-1968B-22G
VB-1968C-23G



SIMCOE
Irrigating-aspirating
cannula with silicone
tubing 23 gauge
VB-1974 reverse

Scanned with OKEN Scanner

CANNULAES HANDLE



KNOLLE
Anterior chamber irrigating cannula 2mm angled tip front opening
VB-2022 (A) 30gauge
VB-2022 (B) 23gauge



GIMBEL
Hydro-dissection Cannula smooth blunt tip front opening 30gauge
VB-2028



AIR INJECTION
Cannula angled 5mm from tip
VB-2034 (A) 26gauge
VB-2034 (B) 27gauge
VB-2034 (C) 30gauge



TROUTMAN
Alpha-chyrnotrypsin cannula delicate olive tip 25gauge
VB-2040



KUGLEN
Irrigating lens manipulator disc shaped tip with 0.25mm diameter peg protruding from its center, 25gauge
VB-2046 (A) Straight
VB-2046 (B) Angled



YOUENS
Lens manipulating cannula front opening angled 30gauge
VB-2052



WELSH
Flat olive tip cannula 25gauge
VB-2058 (A) Straight
VB-2058 (B) Right
VB-2058 (C) Left



DREWS
Posterior capsule polisher slightly angled, 2mm i.d. loop with sandblasted base 0.2 irrigating port facing into circle 25gauge
VB-2064



JENSEN
Posterior capsule polisher carbide impregnated on posterior surface of olive tip front opening gently curved
VB-2070 (A) 25gauge
VB-2070 (B) 23gauge



JENSEN
Capsule polisher carbide impregnated all around olive tips front opening gently curved
VB-2076 (A) 25gauge
VB-2076 (B) 23gauge



KNOLLE
Posterior capsule polisher 2mm angled tip, front opening finely sandblasted
VB-2082 (A) 23gauge
VB-2082 (B) 30gauge



KRATZ
Scratcher front opening with carbide impregnated tips slightly angled 22gauge
VB-2088



SHAPIRO
Posterior capsule polisher carbide impregnated tip with dual 0.3mm side ports curved 23gauge
VB-2094

Scanned with OKEN Scanner

CANNULAES HANDLE



SIMCOE

Cortex extractor cannula
smooth closed tip, with side
opening facing info curved
23 gauge sandblasted
VB-2100 (A) 13mm
VB-2100 (B) 21mm



KELMAN

Cystotome with
sharp point and
sharp cutting edge 22gauge
VB-2136



SIMCOE

Cortex extractor and capsule
polisher front opening facing
into curved 23gauge sandblasted
VB-2106 (A) 13mm
VB-2106 (B) 21mm



KNOLLE-KELMAN

Cystotome small blade with
sharp point and sharp cutting
edge 22gauge
VB-2142



ANIS

Cortex aspirating cannula with
0.4mm side port finely
sandblasted tip 23gauge
VB-2112 (A) Straight
VB-2112 (B) Angled



Mc INTYRE

Cystotome angled 10mm from
tip with sharp point, 23gauge
VB-2148



GRAETHER

Collar button Iris retracting &
capsular polisher cannula
VB-2118 (A) Straight
VB-2118 (B) Angled




SILICONE BULB

With adaptor
VB-2154



DALJIT

Blast cannula
VB-2124



TEPLON HANDLE

VB-2160



HESLIN

Infusion cannula with silicone
tubing and adaptor
VB-2130 (A) 23gauge

Scanned with OKEN Scanner

CAUTERYS, LOOPS & SPATULAS



LENS EXPRESSOR
Hook & Loop
VB-0468



KIRBY
Lens expressor
hook & spatula
VB-0474



WILLS
Hospital cautery
cooper ball
VB-0480



WILLS
Cautery ball point ss
VB-0486



JOSHI
Cautery copper Tip
VB-0492



SNELLEN
Lens loop Serrated
VB-0504

Scanned with OKEN Scanner

CAUTERYS, LOOPS & SPATULAS



WEBER
Lens loop
VB-0510 (A) Small
VB-0510 (B) Large



DASTOOR
Pupil & Iris repositor
VB-0516 (A) Signle
VB-0516 (B) Double



JAEGER
Lid plate
VB-0522



SMITH-FISHER
Iris spatula double
VB-0528



WECKER
Iris spatula
VB-0534



GREEN
Spatula double end
VB-0540

Scanned with OKEN Scanner

CHOPPERS





PHACO CHOPPER
Universal
VB-1386






PHACO CHOPPER
Ball tip
VB-1392






PHACO CHOPPER
Nagahara
VB-1398





COMBO CHOPPER
Mohanraj
VB-1400



PRECHOPPER WITH
Substainer
VB-1404 (A) Straight
VB-1404 (B) Angled
VB-1404 (C) Straight (hole)
VB-1404 (D) Angled (hole)

Scanned with OKEN Scanner

CLAMPS



WORTH
Advancement forcep
VB-0936 (A) Right
VB-0936 (B) Left



SNELLEN
Entropium forcep
VB-0942 (A) Right Big
VB-0942 (B) Left Big

BERKE
Ptosis clamp
VB-0948 20mm



AYER
Chalazion forcep
VB-0954



LAMBERT
Chalazion forcep
VB-0960

Scanned with OKEN Scanner

CLAMPS



BAIRD
Chalazion forcep
VB-0966 8x12mm



HEATH
Chalazion forcep
VB-0972 12x14mm



FRANCIS
Chalazion forcep
VB-0978 12x14mm oval



DESMARRES
Chalazion forcep
VB-0984 20mm



SWISS
Model muscle forcep
VB-0990

Scanned with OKEN Scanner

CURETTES, PROBES, SPOONS & DILATORS




MEYHOEFER
Chalazion curette sharp
VB-0246 (A) 1mm
VB-0246 (B) 2mm
VB-0246 (C) 3mm



HARMS
Trabeculotomy probe
VB-0252 Right & Left



WILDER
Lacrimal dilator
VB-0258



BOWMAN
Lacrimal Probe set
VB-0264 (4 pes set)



MULE
Evisceration spoon
VB-0270



WELLS
Enucleation spoon
VB-0276

Scanned with OKEN Scanner

CURETTES, PROBES, SPOONS & DILATORS



BUNGE
Evisceration spoon
VB-0282



GRAEFE
Lens spoon
VB-0288



WEST
Bone Chisel
VB-0294



WEST
Bone gouge
VB-0300



FREER
Twin elevator sharp & blunt
VB-0306

Scanned with OKEN Scanner

DISSECTORS & CLAMPS



TRAQUIAR'S
Periostoral elevator
double-end
VB-0312



DASTOOR
Lacrimal sac dissector
VB-0318



LANG
Lacrimal sac dissector
VB-0324



PIGTAIL PROBE
With suture hole
VB-0330



SERREFINE
VB-0336 (A) Straight
VB-0336 (B) Curved



BACKHAUS
Towel clamp
VB-0342

TOWEL CLIP
Cross action
VB-0348


Scanned with OKEN Scanner

FORCEPS






EYE DRESSING
Forcep serrated
VB-0618 (A) Staight
VB-0618 (B) Curved





IRIS FORCEP
1x2 teeth
VB-0624 (A) Staight
VB-0624 (B) Curved





SWISS
Model iris forcep
curved 1x2 teeth
VB-0630






HESS
Iris forcep
VB-0636 Swiss Model






BONN
Iris forcep
VB-0642 with platform





BONN
Iris forcep
VB-0648 with platform

Scanned with OKEN Scanner





BISHOP-HARMON
Iris forcep
VB-0654 (A) 1x2 teeth
VB-0654 (B) serrated





KIRBY
Iris forcep 1x2 teeth
VB-0660





UTILITY
Forcep
VB-0666 1x2 straight





BONACCOL-TO
Utility forcep
VB-0672





WILLS HOSPITAL
Utility forcep
VB-0678





BARRAQUER
Colibri forcep
VB-0684 0.2mm teeth

Scanned with OKEN Scanner

FORCEPS





CASTROVIEJO
Colibri forcep
very delicate 1x2 teeth
VB-0690





LIM'S
Corneal forcep 1x2 teeth
VB-0708 (A) Small
VB-0708 (B) Big





CASTROVIEJO
Suture forcep teeth
VB-0720 0.5mm Straight





CASTROVIEJO
Wide handle 1x2
teeth forcep
VB-0726





PAUFIQUE
Suture forcep 1x2 teeth
VB-0732



ALABAMA
Suture forcep
VB-0738

Scanned with OKEN Scanner


FORCEPS



JAYLE'S
Suture forcep 1x2 teeth
VB-0744



MOORFIELD
Suture and conjunctival forcep
VB-0750





FIXATION
Forcep
VB-0756 (A) 1x2 Teeth
VB-0756 (B) 2x3 Teeth
VB-0756 (C) 3x4 Teeth



GRAEFE
Fixation forcep
VB-0762





DASTOOR
Superior-rectus forcep
VB-0792






TROUTMAN
Superior-rectus forcep
VB-0798

Scanned with OKEN Scanner

FORCEPS






PIERSE TYPE
Micro forcep
VB-0804 0.2mm tip





PIERSE TYPE
Micro forcep curved
VB-0810 0.2mm tip





PIERSE TYPE
Fine forcep
VB-0816 0.2mm tip





PIERSE
Colibri Micro forcep
VB-0822 0.2mm tip



ELSCHNIG
Capsular forcep
VB-0840



BARRAQUER
Cilia forcep slightly angled
VB-0852

Scanned with OKEN Scanner

FORCEPS



ZIEGLER
Cilia forcep
VB-0858





BARRAQUER
Tying forcep
VB-0864 (A) Straight
VB-0864 (B) Curved





HARMS
Tying forcep smooth jaws
VB-0870 (A) Straight
VB-0870 (B) Curved





SINSKEY
Tying forcep curved
VB-0876





TENNANT
Tying forcep extra delicate jaws
VB-0882 (A) Straight
VB-0882 (B) Curved





MICRO
Micro tying forcep
VB-0888 (A) Straight Small
VB-0888 (B) Curved Small

Scanned with OKEN Scanner

FORCEPS





MICRO
Micro tying forcep
VB-0894 (A) Straight Big
VB-0894 (B) Curved Big



CASTROVIEJO
Wide handle tying forcep
VB-0900





BONN
Micro tying forcep
VB-0906 Short





JAYLE'S
Tying forcep
VB-0912 Straight



JEWELERS
Forcep
VB-0918





JEWELERS
Forcep angled
VB-0924 (A) Straight
VB-0924 (B) Curved

Scanned with OKEN Scanner

HOOKS & POLISHERS





HIRSCHMAN
Iris push & pull hook
VB-1320 (A) Straight
VB-1320 (B) Angled





NIGHTINGALE
Capsule polishing
curette 1.7mm
VB-1338





NIGHTINGALE
Capsule polishing Curette
VB-1344 (A) Right
VB-1344 (B) Left



NORDAN
Capsule elevator
for in-th-bag
IOL placement
VB-1350





KANSAS
Nucleus trisector
VB-1356 (A) Right
VB-1356 (B) Left





KANSAS
Nucleus bisector
VB-1362

Scanned with OKEN Scanner

HOOKS



ARRUGA
Hook
VB-0402




CULLER
Muscle hook
VB-0408



GREEN
Strabismus hook
VB-0414



JAMESON
Strabismus hook
VB-0420 (A) Small
VB-0420 (B) Big



GRAEFE
Strabismus hook
VB-0426

Scanned with OKEN Scanner

HOOKS



O'CONNOR
Tenotomy hook
VB-0432



GASS
Detachment hook
With oval hole
VB-0438



FLAT
Hook
VB-0444



CHAVASSE'S
Muscle hook
VB-0450



SMITH
Lens expressor
VB-0456



KIRBY
Lens expressor
Hook double
VB-0462

Scanned with OKEN Scanner

IOL MANIPULATORS



BECHERT
Lens nucleus rotator
blunt tip angled
VB-1206 (A) 7mm from tip



CLAYMAN
Nucleus rotator
v-notched tip
VB-1212



JAFFE
Lens spatula
flat smooth
VB-1230



GRAETHER
Collar button
micro iris retractor
VB-1236 (A) Straight
VB-1236 (B) Angled



JAFFE-KNOLLE
Iris hook very blunt
small angled
VB-1242



SINSKEY
Lens manipulating hook
VB-1248 Angled



LEWICKY
Lens manipulating hook
angled
VB-1254

Scanned with OKEN Scanner

IOL MANIPULATORS





DRYSDALE
Nucleus manipulator
paddle shaped tip
VB-1260





KUGLEN
Iris hook and lens manipulator
push-pull model
VB-1266 Angled





FENZL
Lens manipulating hook
T shaped tip
VB-1272 Angled





KNOLLE
Iris hook very blunt
VB-1284 Straight





BONN
Micro iris hook
VB-1290





PISACANO
Nucleus rotator
flat blunt tip
VB-1308





MANN-HUNKELER
Capsule retractor
with 0.3mm diameter
semi-ball
VB-1314

Scanned with OKEN Scanner

KERATOMES KNIVES SPUDS



CASTROVIEJO
Keratome
VB-0168



ZIEGLER
Iris knife
VB-0174



TOOK KNIFE
Corneal knife
VB-0180



GILLS
Corneal knife
VB-0186



PAUFIQUE
Keratoplasty Knife
VB-0192



BARRAQUER
Keratoplasty Knife
VB-0198

Scanned with OKEN Scanner

KERATOMES KNIVES & SPUDS



GRAEFE
Cystotme
VB-0216



GRAEFE
Cystotme & spoon
VB-0222



DIX
Foreign body spud
VB-0228




ROLLET
Rugine for lacrimal sac
VB-0234



HOCKEY KNIFE
Corneal Dissector
VB-0240

Scanned with OKEN Scanner

MALLETS & FORCEPS



MALLET
For DCR
VB-0354



HARTMANN
Mosquito forcep 3.5”
delicate
VB-0360 (A) Straight
VB-0360 (B) Curved



HALSTED
Mosquito hemostatic
forceps 5.5"
VB-0366 (A) Straight
VB-0366 (B) Curved

Scanned with OKEN Scanner

MARKERS FIXATION RING





HOFFER
Optic zone marker
With cross hairs
VB-0582 (A) 5.0 mm
VB-0582 (B) 5.5 mm
VB-0582 (C) 6.0 mm



THORNTON
Fixation ring with
Atraumatic teeth 16 mm
VB-0606



HOFMAN-THORNTON
Global fixation ring with
VB-0612 (A) 16mm
VB-0612 (B) 13mm



MANDEZ
Degree Gauge Marker
VB-0614

Scanned with OKEN Scanner

PUNCHES & TREPHINES



KERRISON

D C R Punch

VB-0372 (A) 1.0mm

VB-0372 (B) 2.0mm

VB-0372 (C) 3.0mm

VB-0372 (D) 4.0mm



KELLY'S

Detachment Punch

VB-0378



GASS

Trabeculotomy punch

VB-0384



CASTROVIEJO

Trephine

6mm to 11mm

VB-0390A - Disposable

VB-0390B - Stainless Steel

Scanned with OKEN Scanner

RHEXIS, HOLDER & FOLDER



Mc PHERSON
Forcep
VB-1002 (A) 5mm
VB-1002 (B) 8mm



Mc PHERSON
Forcep delicate 11mm angled
VB-1014



KELMAN-Mc PHERSON
Forcep
VB-1032 (A) 5mm
VB-1032 (B) 8mm



KELMAN-Mc PHERSON
Forceps 11mm plain
Angled delicate
VB-1038





KRAFF
Capsule tag forcep
11mm long jaws
VB-1044





UTRATA
Capsularhexis forcep
delicate grasping tips
VB-1050

Scanned with OKEN Scanner

RHEXIS, HOLDER & FOLDER





MASKET
Capsularhexis forcep
VB-1056



SCLERAL
Plug Holding Forceps
VB-1068





DODICK
Nucleus cracker
cross-action
VB-1074



KANSAS
Nucleus removal forcep
VB-1080





DALJIT
Lens holding forcep
VB-1086





ARRUGA
Capsular forcep
VB-1104 Standard

Scanned with OKEN Scanner

RHEXIS, HOLDER & FOLDER



MICRO
Ring forcep
VB-1110



BLAYDES
Lens holding forcep
VB-1116




ANIS
Lens holding forcep
extra delicate
VB-1122





KNOLLE-SHEPARD
Lens holding forcep
with peg serrated jaws
VB-1128





SHEPARD
Lens holding forcep
with serrated lower jaws
VB-1134 without lock



IRIS CLAW
Lens holding forcep
VB-1140

Scanned with OKEN Scanner

RETRACTORS & KNIVES



DESMARRES
Lid retractor
VB-0108



PAUL
Retractor, blunt
VB-0114



CATPAW
Retractor 4
prongs for sac
VB-0120



SCHEPENS
Orbital retractor
VB-0126



DASTOOR
Iris retractor
VB-0132



GRAEFE
Cataract knife
VB-0150

Scanned with OKEN Scanner

SCISSORS



ENUCLEATION

Scissor

VB-1524 (A) Half Curved

VB-1524 (B) Fll Curved

VB-1524 (C) St blunt tip



IRIS

Scissor 9 cm

VB-1530 (A) Straight

VB-1530 (B) Curved



IRIS

Scissor 10.5 cm

VB-1536 (A) Straight

VB-1536 (B) Curved



STEVENS

Tenotomy scissor

VB-1548 (A) Straight

VB-1548 (B) Curved

Scanned with OKEN Scanner

SCISSORS



NOYES
Iris scissor
VB-1554 (A) Blunt
VB-1554 (B) Sharp



VANNAS
Scissor sharp tips
VB-1560 (A) Straight blade 7mm
VB-1560 (B) Straight blade 10mm



VANNAS
Scissor sharp tips
VB-1566 (A) Curved blade 7mm
VB-1566 (B) Curved blade 10mm



VANNAS
Scissor sharp tips
VB-1572 angled to side



CAPSULOTOMY
Scissor curved blunt tips
VB-1578 14mm blade



GILLS – VANNAS
Capsulotomy scissor sharp tips
VB-1584
Angled forward & 11mm

Scanned with OKEN Scanner

SCISSORS



STERN-GILLS
Capsulotomy scissor sharp tips
VB-1590 Curved



STERN-GILLS
Capsulotomy scissor
pointed blade
VB-1596 angled 45° at 10mm tips



WESTCOTT
Tenotomy scissor curved
blunt tips standard blade
VB-1602



WESTCOTT
Tenotomy scissor curved
blunt tips medium blade
VB-1608



WESTCOTT
Conjunctival scissor curved
blunt tips small blade
VB-1614



WESTCOTT
Stitch scissor
standard blade
VB-1620

Scanned with OKEN Scanner


SCISSORS



WESTCOTT
Stitch scissor
sharp pointed tips
medium blade
VB-1626



CASTROVIEJO
Corneal scissor
angled to side blunt tips
VB-1632 small blade 11mm



IRIS
Scissor sharp tips
small blade
VB-1638 (A) Straight
VB-1638 (B) Curved



MICRO
Corneal Scissor
blunt tips small blade
VB-1650 Half Curved



MICRO
Corneal scissor
blunt tips small blade
VB-1656 Full Curved



MICRO
Corneal scissor
blunt tips small blade (hole)
VB-1668 Slightly Curved

Scanned with OKEN Scanner

SCISSORS



MICRO
Corneal scissor
blunt tips small blade (hole)
VB-1674 Half Curved



MICRO
Corneal scissor
blunt tips small blade (hole)
VB-1680 Full Curved



CASTROVIEJO
Corneal scissor
blunt tips large blade
VB-1686 Slight Curved



CASTROVIEJO
Corneal scissor
blunt tips large blade
VB-1692 Half Curved




CASTROVIEJO
Corneal scissor
blunt tips large blade
VB-1698 Full Curved



CASTROVIEJO
Corneal scissor
VB-1710 (A) Right Medium
VB-1710 (B) Left Medium

Scanned with OKEN Scanner


SPECULUMS




BARRAQUER
Wire speculum
VB-0006 (A) Adult
VB-0006 (B) Child



BARRAQUER
Wire eye speculum
VB-0012



KRATZ-BARRAQUER
Wire speculum
VB-0024



LIEBERMAN
Eye speculum
VB-0030 (A)
Phaco Speculum
VB-0030 (B) Solid Blade




CLARK
Eye speculum
VB-0036

Scanned with OKEN Scanner

SPECULUMS




LANCASTER
Eye speculum
VB-0042



MURDOCK
Eye speculum
VB-0048




SAUER
Eye speculum
VB-0054



INFANT
Wire speculum
VB-0060



ALPHONSO
ROP Speculum
VB-0064



WEISS
Eye speculum
VB-0066

Scanned with OKEN Scanner

SPECULUMS




SCLERAL DEPRESSOR
VB-0078 A - Single End
VB-0078 B - Double End



SCHEPENS
Scleral depressor
VB-0084 (A) 8mm
VB-0084 (B) 10mm



MUELLER'S
Lacrimal Sac Retractor
VB-0090



STEVENSON
Lacrimal sac retractor
VB-0096



JAFFE
Lid retractor
VB-0102 (A) Small
VB-0102 (B) Large

Scanned with OKEN Scanner

SPATULAS, CALIPERS & MARKERS



SPATULA
With needle
VB-0546



CASTROVIEJO
Cyclodialysis spatula
VB-0558 (A) Double
VB-0558 (B) Elsching Single



CASTROVIEJO
Caliper 20 mm
VB-0564 (A) Straight
VB-0564 (B) Curved



OSHER-NEUMANN
Radial marker
VB-0576 (A) 4 Blades
VB-0576 (B) 6 Blades
VB-0576 (C) 8 Blades
VB-0576 (D) 12 Blades
VB-0576 (E) 16 Blades

Scanned with OKEN Scanner

STERILIZATION CASE



INFUSION
Handle
VB-2166 (A) 9.0 long
VB-2166 (B) 12.5 long



TREPHINE PUNCH
VB-2172



S S BOX
Metal tray with one silicone mat
VB-2178

Metal tray with two silicone mat
VB-2184 (A) Height
VB-2184 (B) Width

Metal tray with four silicone mat
VB-2190

Scanned with OKEN Scanner

VITREO - RETINAL INSTRUMENTS & PEELARS

DORC TYPE FORCEPS & SCISSOR





VB-1410
Straight forceps



VB-1416
Serrated Forceps




VB-1422
End gripping



VB-1428
45° Forceps



VB-1430
IOL Holding
Forceps 20G



VB-1432
ILM Forcep



VB-1434
Forceps



VB-1440
Straight Scissor
Sharp tip



VB-1446
Curved Scissor
Sharptip



VB-1448
IOL Cutting
Scissor 20G



VB-1452
Vertical Scissor
45°



VB-1464
Capsulorhexis 45°



VB-1470
Capsulorhexis 90°



VB-1476
Foreign body

MEMBRANE PEELAR





VB-1482
Peelar



VB-1488
Spatula curved
pointed tip



VB-1494
Scratcher clow
type



VB-1500
Spatula hooked



VB-1506
Spatula hooked
short length



VB-1512
Scratcher



VB-1518
Spatula knob

Scanned with OKEN Scanner

VIAN

Vision beyond limits

Vian Biotech Private Limited

534, Water Tank Road,
Rajarajeshwari Nagar, Bangalore - 560 098.
Mob.: +91 97405 26706 / 99168 60403
Email : info@vianbiotech.com
www.vianbiotech.com

Scanned with OKEN Scanner